



Guidance for the selection of Geister SurgiBOX Micro Type

Type 1

Outer dimensions:	220 x 165 mm
Inner dimensions:	200 x 140 mm
Maximum instrument length:	200 mm



77-0101

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 3 pcs
- Forceps (Standard and Micro) max. 5 pcs
- Nerve hooks max. 7 pcs



77-0111

Options for content e.g.:

- Same dimensions and content options than 77-0101 **with** twice the content.



77-0121

Options for content e.g.:

- Silicone benches can be customized.

Type 2

Outer dimensions:	420 x 125 mm
Inner dimensions:	400 x 100 mm
Maximum instrument length:	100 mm



77-0102

Options for content e.g.:

- Microinstruments for ophthalmology max. 15 pcs



77-0112

Options for content e.g.:

- Same dimensions and content options than 77-0102 **with** twice the content



77-0122

max. instrument length: 400 mm

Options for content e.g.:

- Rigid / semiflexible optics max. 2 pcs
- ValveGate max. 2-3 pcs.

Type 3

Outer dimensions:	420 x 175 mm
Inner dimensions:	400 x 150 mm
Maximum instrument length:	150 mm
(- 165 mm with more diagonal angle)	



77-0103

Options for content e.g.:


- Microinstruments max. 15 pcs



77-0113

Options for content e.g.:

- Same dimensions and content options than 77-0103 **with** twice the content.



77-0123

Options for content e.g.:

- Silicone benches can be customized.

77-0103V

Options for content e.g.:

- Rigid / semiflexible optics max. 3 pcs
- ValveGate max. 5 pcs.

77-0113V

Options for content e.g.:

- Same dimensions and content options than 77-0103V **with** twice the content.

Type 4

Outer dimensions:	300 x 220 mm
Inner dimensions:	280 x 195 mm
Maximum instrument length:	195 mm



77-0104

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 6 pcs
- Forceps (Standard and Micro) max. 10 pcs
- Nerve hooks max. 14 pcs



77-0114

Options for content e.g.:

- Same dimensions and content options than 77-0104 **with** twice the content.



77-0124

Options for content e.g.:

- Silicone benches can be customized.

Type 5

Outer dimensions:	285 x 280 mm
Inner dimensions:	265 x 255 mm
Maximum instrument length:	255 mm



77-0105

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 6 pcs
- Forceps (Standard and Micro) max. 10 pcs
- Nerve hooks max. 14 pcs
- Ring instruments max. 5 pcs



77-0115

Options for content e.g.:

- Same dimensions and content options than 77-0105 **with** twice the content.




77-0125

Options for content e.g.:

- Silicone benches can be customized (e.g. for small retractor).

Type 6

Outer dimensions:	420 x 125 mm
Inner dimensions:	400 x 100 mm
Maximum instrument length:	100 mm



77-0106

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 9 pcs
- Forceps (Standard and Micro) max. 12-14 pcs
- Nerve hooks max. 20 pcs
- Ring instruments max. 7 pcs



77-0116

Options for content e.g.:

- Same dimensions and content options than 77-0106 **with** twice the content.



77-0126

Options for content e.g.:

- Silicone benches can be customized (e.g. for standard instruments or small retractors).

Type 7

Outer dimensions:	400 x 290 mm
Inner dimensions:	380 x 265 mm
Maximum instrument length:	265 mm (- 280 mm with more diagonal angle)



77-0107

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 9 pcs
- Forceps (Standard and Micro) max. 12-14 pcs
- Nerve hooks max. 20 pcs
- Ring instruments max. 7 pcs



77-0117

Options for content e.g.:

- Same dimensions and content options than 77-0107 **with** twice the content.



77-0127

Options for content e.g.:

- Silicone benches can be customized (e.g. for standard instruments or small retractors).

Type 8

Outer dimensions:	450 x 260 mm
Inner dimensions:	430 x 230 mm
Maximum instrument length:	230 mm



Example

77-0108

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 10 pcs
- Forceps (Standard and Micro) max. 15 pcs
- Nerve hooks max. 25 pcs
- Ring instruments max. 8 pcs



77-0118

Options for content e.g.:

- Same dimensions and content options than 77-0108 **with** twice the content.




77-0128

Options for content e.g.:

- Silicone benches can be customized (e.g. for standard instruments or small retractors).

Type 9


Outer dimensions:	400 x 150 mm
Inner dimensions:	380 x 125 mm
Maximum instrument length:	125 mm



77-0109

Options for content e.g.:

- Microinstruments for ophthalmology max. 14 pcs



77-0119

Options for content e.g.:

- Same dimensions and content options than 77-0109 **with** twice the content.



77-0129

Options for content e.g.:

- Silicone benches can be customized.

77-0109V

Options for content e.g.:

- Rigid / semiflexible optics max. 2-3 pcs
- ValveGate max. 3 pcs.

77-0119V

Options for content e.g.:

- Same dimensions and content options than 77-0109V **with** twice the content.

Type 10

Outer dimensions:	495 x 270 mm
Inner dimensions:	470 x 250 mm
Maximum instrument length:	250 mm



77-0110

Options for content e.g.:

- Microinstruments (Needle driver, scissor) max. 12 pcs
- Forceps (Standard and Micro) max. 18 pcs
- Nerve hooks max. 30 pcs
- Ring instruments max. 9 pcs



77-0120

Options for content e.g.:

- Same dimensions and content options than 77-0110 **with** twice the content.



77-0130

Options for content e.g.:

- Silicone benches can be customized (e.g. for standard instruments or small retractors).